

LISTEN.
THINK.
SOLVE.®

PharmaSuite®



NON-FUNCTIONAL REQUIREMENTS

RELEASE 8.4

FUNCTIONAL REQUIREMENT SPECIFICATION

PUBLICATION PSFRSNF-RM001E-EN-E-DECEMBER-2017

Supersedes publication PSFDNF-RM001D-EN-E



Allen-Bradley • Rockwell Software

Rockwell
Automation

Contact Rockwell See contact information provided in your maintenance contract.

Copyright Notice © 2017 Rockwell Automation Technologies, Inc. All rights reserved.
This document and any accompanying Rockwell Software products are copyrighted by Rockwell Automation Technologies, Inc. Any reproduction and/or distribution without prior written consent from Rockwell Automation Technologies, Inc. is strictly prohibited. Please refer to the license agreement for details.

Trademark Notices FactoryTalk, PharmaSuite, Rockwell Automation, Rockwell Software, and the Rockwell Software logo are registered trademarks of Rockwell Automation, Inc.

The following logos and products are trademarks of Rockwell Automation, Inc.:

FactoryTalk Shop Operations Server, FactoryTalk ProductionCentre, FactoryTalk Administration Console, FactoryTalk Automation Platform, and FactoryTalk Security.
Operational Data Store, ODS, Plant Operations, Process Designer, Shop Operations, Rockwell Software CPGSuite, and Rockwell Software AutoSuite.

Other Trademarks ActiveX, Microsoft, Microsoft Access, SQL Server, Visual Basic, Visual C++, Visual SourceSafe, Windows, Windows 7 Professional, Windows Server 2008, Windows Server 2012, and Windows Server 2016 are either registered trademarks or trademarks of Microsoft Corporation in the United States and/or other countries.

Adobe, Acrobat, and Reader are either registered trademarks or trademarks of Adobe Systems Incorporated in the United States and/or other countries.

ControlNet is a registered trademark of ControlNet International.

DeviceNet is a trademark of the Open DeviceNet Vendor Association, Inc. (ODVA).

Ethernet is a registered trademark of Digital Equipment Corporation, Intel, and Xerox Corporation.

OLE for Process Control (OPC) is a registered trademark of the OPC Foundation.

Oracle, SQL*Net, and SQL*Plus are registered trademarks of Oracle Corporation.

All other trademarks are the property of their respective holders and are hereby acknowledged.

Warranty This product is warranted in accordance with the product license. The product's performance may be affected by system configuration, the application being performed, operator control, maintenance, and other related factors. Rockwell Automation is not responsible for these intervening factors. The instructions in this document do not cover all the details or variations in the equipment, procedure, or process described, nor do they provide directions for meeting every possible contingency during installation, operation, or maintenance. This product's implementation may vary among users.

This document is current as of the time of release of the product; however, the accompanying software may have changed since the release. Rockwell Automation, Inc. reserves the right to change any information contained in this document or the software at any time without prior notice. It is your responsibility to obtain the most current information available from Rockwell when installing or using this product.

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Chapter 1	Introduction	1
	Typographical Conventions	2
Chapter 2	Performance Requirements (SR1095.10+)	3
	Master Recipe - Performance (SR1095.10.4+)	3
	Performance of Master Recipe Import 1 (SR1095.10.4.1)	4
	Performance of Master Recipe Export 1 (SR1095.10.4.2)	4
	Performance of Master Recipe Saving 1 (SR1095.10.4.3)	4
	Performance of Prepare for Status Change 1 (SR1095.10.4.6)	4
	Performance of Master Recipe Report 1 (SR1095.10.4.5)	4
	Performance of Order Explosion 1 (SR1095.10.4.4)	4
	Change Request - Performance (SR1095.10.1+)	5
	Performance of Change Request Execution 1 (SR1095.10.1.1)	5
	Performance of Change Request Execution 2 (SR1095.10.1.2)	6
	Performance of Usage List Creation 1 (SR1095.10.1.3)	6
	Equipment - Performance (SR1095.10.2+)	6
	Performance of Unit Class Export 1 (SR1095.10.2.1)	7
	Performance of Unit Class Import 1 (SR1095.10.2.2)	7
	Performance of Filter Class Export 1 (SR1095.10.2.5)	7
	Performance of Filter Class Import 1 (SR1095.10.2.6)	8
	Recipe Execution - Performance (SR1095.10.7+)	8
	Performance of Recipe Execution 1 (SR1095.10.7.1)	8
	Printing of Barcode Label and Dispensing Report - Performance	9
	Performance of Printing (SR1090.4.6)	9
	EBR Server - Performance (SR1095.10.3+)	9

Performance of Order Processing 1 (SR1095.10.3.1)	9
TOM Server - Performance (SR1095.10.5+).....	10
Performance of ETO Processing 1 (SR1095.10.5.1)	10
Exception Dashboard - Performance (SR1095.10.8+)	10
Performance of Review by Exception 1 (SR1095.10.8.1).....	11
Performance of Batch Report Printing 1 (SR1095.10.8.2).....	11
Chapter 3 Compliance Requirements (SR1095.50+).....	13
Audit Trail (SR1095.50.1)	13
Electronic Signatures (SR1095.50.2).....	13
Single and Double Electronic Signatures (SR1095.50.2.1).....	14
Configuring Electronic Signatures (SR1095.50.2.2)	15
Performing Electronic Signatures (SR1095.50.2.3).....	15
Data of an Electronic Signature (SR1095.50.2.4).....	16
Access Control (SR1095.50.3).....	16
Authentication for System Login and Electronic Signature (SR1095.50.3.1)	16
Management of Users (SR1095.50.3.2)	18
Management of User Groups (SR1095.50.3.3).....	20
Management of Access Privileges (SR1095.50.3.4).....	20
Password Security (SR1095.50.3.5)	20
Password Change (SR1095.50.3.6)	21
Chapter 4 Internationalization and Localization (SR1095.1+)	23
Skills for Localization (SR1095.1.1).....	23
Tool Support (SR1095.2.4)	23
Internationalized Output Media (SR1095.2)	23
Locale-related Requirements.....	24
Supported Locales (SR1095.1.2).....	24
Provided Locales (SR1095.2.3+)	24
Locales within One Installation (SR1095.2.1)	24
Choice of Locale (SR1095.2.2).....	25
Configurable Locale-specific Settings (SR1095.1.3).....	25

Chapter 5	Installation, Backup, and Upgrade (SR1095.4+)	27
	Installation of PharmaSuite (SR1095.4.1).....	27
	Initial Installation (SR1095.4.1.1).....	27
	Existing Database (SR1095.4.1.7)	27
	Installation Result (SR1095.4.1.8)	27
	Scope of Express Edition Installer (SR1095.4.1.6).....	28
	Installation Manual (SR1095.4.1.4)	28
	Backup (SR1075.6).....	28
	Backup of a Running System (SR1075.6.1).....	28
	Backup and Restoring (SR1075.6.3).....	28
	Upgrade Tool (SR1075.8+)	28
	Basic Upgrade Engine (SR1075.8.1).....	29
	System Update Tasks (SR1075.8.2+).....	29
	System Update and Data Migration Packages (SR1075.8.3+)	31
	Extensibility (SR1075.8.5+)	32
	Pre-install Configuration (SR3146.7)	32
Chapter 6	User Experience (SR1095.5+).....	33
	Messages (SR1095.5.1.1)	33
	Device-specific Usability (SR1095.5.1.4+)	33
	Touch Screen (SR1095.5.1.4.2).....	33
	Virtual Keyboard (SR1095.5.1.4.2.1)	33
	Look and Feel (SR1095.5.7).....	33
	GUI Design (SR1095.5.7.3)	33
	GUI Adaptation (SR1095.5.7.1).....	33
	Dialog Boxes (SR1095.5.7.2)	33
	Grid Components (SR1095.5.7.4)	34
	Extensibility (SR1095.5.1.2+).....	34
	System - Customization (SR1095.5.1.2.1)	34
	System - View Model-related Objects (SR1095.5.1.2.3).....	34
	System - Material Attribute (SR1095.5.1.2.9).....	34

System - ERP BOM Attribute (SR1095.5.1.2.11).....	34
System - Batch Status (SR1095.5.1.2.10)	35
Master Recipe - Approval Workflow (SR1095.5.1.2.13).....	35
Expression Editor - Context Data (SR1095.5.1.2.14).....	35
Building Block - Name Spaces (SR1095.5.1.2.4)	35
Building Block - Mainstream GUI Programming (SR1095.5.1.2.6)	35
Execution - Autostart (SR1095.5.1.2.2)	35
Execution - Select Order Grid (SR1095.5.1.2.5)	36
Execution - EBR Header Bar (SR1095.5.1.2.8)	36
Reporting (SR1095.5.1.2.16+)	36
Precision of Numerical Data (SR1095.5.1.3+)	36
Fixed Point Numbers (SR1095.5.1.3.1).....	36
Floating Point Arithmetic (SR1095.5.1.3.2)	37
Errors at Calculation (SR1095.5.1.3.3).....	37
Input Data Errors (SR1095.5.1.3.4)	37
Rounding (SR1095.5.1.3.5).....	37
Concurrent Data Access - Pessimistic Locking (SR1095.5.1.5)	37
Restart and Recovery (SR1095.8)	37
Chapter 7 Reference Documents	39
Chapter 8 Document Information	41
Approval	41
Version Information	41
Revision History.....	41
Index	47

Figure 1: Example of a single electronic signature	14
Figure 2: Example of a double electronic signature to support witness role	14

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Introduction

This document details the non-functional requirements of PharmaSuite.

Related requirements that already existed before PharmaSuite 8.4 have been moved from the following document to this Functional Requirement Specification:

- Detailed Requirements Specification (DRS) document "PharmaSuite 8.0 - DRS Information Management" (106810)

Each requirement is composed of a name and a unique identifier (e.g. Performance of Master Recipe Import 1 (SR1095.10.4.1)). If a requirement's meaning is for requirement grouping only, the identifier is appended by a plus sign (e.g. Master Recipe - Performance (SR1095.10.4+)).

In some cases, additional context information is available, indicated in the document by a frame and a gray background color. This context information is related to the respective requirement, but not part of the formal requirement description.

Along with the update of this document, almost all requirement names and descriptions from the DRS documents have been updated slightly. However, this has been an editorial update only, without any impact on the content and scope of the respective requirements. Therefore, those changes are only related to the document within the revision history (page [41](#)).

The revision history (page [41](#)) lists the changes made to the document with PharmaSuite 8.3 as the comparison baseline. Changes related to a requirement are marked as "Editorial", "Update", "New", or "Deleted", changes to the additional context information are marked as "Context information-related".

Typographical Conventions

This documentation uses typographical conventions to enhance the readability of the information it presents. The following kinds of formatting indicate specific information:

Bold typeface	Designates user interface texts, such as
	■ window and dialog titles
	■ menu functions
	■ panel, tab, and button names
	■ box labels
	■ object properties and their values (e.g. status).

Performance Requirements (SR1095.10+)

The detailed list of the PharmaSuite performance requirements is structured in the following sections:

- **Master Recipe - Performance (SR1095.10.4+)** requirements (page [3](#))
- **Change Request - Performance (SR1095.10.1+)** requirements (page [5](#))
- **Equipment - Performance (SR1095.10.2+)** requirements (page [6](#))
- **Recipe Execution - Performance (SR1095.10.7+)** requirements (page [8](#))
- **Printing of Barcode Label and Dispensing Report - Performance** requirements (page [9](#))
- **EBR Server - Performance (SR1095.10.3+)** requirements (page [9](#))
- **TOM Server - Performance (SR1095.10.5+)** requirements (page [10](#))
- **Exception Dashboard - Performance (SR1095.10.8+)** requirements (page [10](#))

For performance requirements of phase building blocks, please refer to

- **Show historical data chart** phase in "Functional Requirement Specification Equipment Automation Phases" [A5] (page [39](#)).
- **Set OPC values** phase in "Functional Requirement Specification Equipment Automation Phases" [A5] (page [39](#)).
- **Get OPC values** phase in "Functional Requirement Specification Equipment Automation Phases" [A5] (page [39](#)).

Master Recipe - Performance (SR1095.10.4+)

For details about recipe management, see "Functional Requirement Specification Recipe and Workflow Management" [A1] (page [39](#)).

Reference scenario 1

The reference recipe holds a total of:

- 3 unit procedures,
- 5 operations,
- 20 phases, and

- 100 process parameters.

Boundary scenario 1

The boundary recipe holds a total of:

- up to 30 unit procedures,
- up to 35 operations,
- up to 3,000 phases, and
- up to 12,000 process parameters.

Its procedural structure includes various branches and loops.

Performance of Master Recipe Import 1 (SR1095.10.4.1)

Based on the Boundary scenario 1 (page 4), the import of the boundary recipe does not take longer than 10 minutes.

Performance of Master Recipe Export 1 (SR1095.10.4.2)

Based on the Boundary scenario 1 (page 4), the export of the boundary recipe does not take longer than 10 minutes.

Performance of Master Recipe Saving 1 (SR1095.10.4.3)

Based on the Boundary scenario 1 (page 4), saving the boundary recipe does not take longer than 10 seconds.

Performance of Prepare for Status Change 1 (SR1095.10.4.6)

Based on the Boundary scenario 1 (page 4), the **Prepare for status change** function does not take longer than 20 seconds.

Performance of Master Recipe Report 1 (SR1095.10.4.5)

Based on the Reference scenario 1 (page 3), the response time for printing the Master Recipe Report of the reference recipe in the **Verification** status (time between "trigger the print function" and "display of the preview") does not take longer than 1 minute.

Performance of Order Explosion 1 (SR1095.10.4.4)

Based on the Boundary scenario 1 (page 4), the order explosion of a boundary recipe-related order does not take longer than 1 minute.

Change Request - Performance (SR1095.10.1+)

For details about change requests, see "Functional Requirement Specification Recipe and Workflow Management" [A1] (page 39).

Reference scenario 1

For the performance test of the change request, the scope of the reference scenario is defined by replacing a **ShowSOP** custom phase within a total of 100 affected reference recipes.

The reference recipe is built with a total of 10 unit procedures, 22 operations, and 192 phases.

The recipe structure is as follows:

- 4 parallel Dispense unit procedures
 - Each with 1 Dispense operation
 - Each with 9 phases
 - Incl. 1 phase based on the **ShowSOP** custom phase
- 6 subsequent processing unit procedures
 - Each with 3 subsequent operations
 - 2 operations
 - Each with 3 subsequent phases
 - Incl. 1 phase based on the **ShowSOP** custom phase
 - 1 operation
 - 20 subsequent phases
 - Incl. 4 phases based on the **ShowSOP** custom phase

Each procedure, unit procedure, and operation of the master recipe is built based on a custom building block in the **Approved** status. In total, 2 custom unit procedures, 4 custom operations, and 14 custom phases are used.

Boundary scenario 1

For the boundary test of the change request, the scope of the boundary scenario is defined by replacing a **ShowSOP** custom phase within a total of 1000 affected reference recipes.

Performance of Change Request Execution 1 (SR1095.10.1.1)

Based on the Reference scenario 1 (page 5), the execution of the change request does not take longer than 40 minutes.

Performance of Change Request Execution 2 (SR1095.10.1.2)

Based on the Boundary scenario 1 (page 5), the execution of the change request is completed without any system errors.

Performance of Usage List Creation 1 (SR1095.10.1.3)

Based on the Boundary scenario 1 (page 5), the initial creation of the usage list within the change request is completed within 1 minute.

The expected response time is based on the following assumptions:

- All recipes that already existed in the database before the reference scenario was created (100 new recipes), are already indexed.
- Along with the creation of the usage list, new indexes are created only for the new recipes of the reference scenario.

Equipment - Performance (SR1095.10.2+)

For details about equipment management, see "Functional Requirement Specification Data Management" [A2] (page 39).

Reference scenario 1

For the performance test of the equipment management functionality, a reference equipment model is considered with a total of 230 classes, 5,000 entities, and 4,225 property types.

The model is structured as follows:

- Units, each with 67 properties
- Filters, hoses, instruments, each with 15 properties
- Number of class-entity-assignments and used property types
(assumption: re-usage of property types between same type of classes of 50 %):
 - 100 unit classes, each with 10 entities
 - Based on 3,250 property types
 - 10 filter classes, each with 100 entities
 - Based on 75 property types
 - 20 hose classes, each with 50 entities
 - Based on 150 property types
 - 100 instrument classes, each with 20 entities
 - Based on 750 property types

Boundary scenario 1

For the boundary test of the equipment management functionality, the reference scenario is scaled up completely by a factor of 10, which means 2,300 classes, 50,000 entities, 42,250 property types.

Update scenario 1

Based on Reference scenario 1 (page 6), an update rate per day of 1 % is assumed, which means that, e.g. 50 entities are updated per day with 1 runtime property being updated per entity.

This results in:

- 18,250 change history entries per year,
- 36,500 change history entries after 2 years, and
- an average of 7-8 update-related change history entries per entity after 2 years.

Update scenario 2

Based on Reference scenario 1 (page 6), an update rate per day of 10 % is assumed, which means that, e.g. 100 entities are updated per day with 5 runtime properties being updated per entity.

This results in:

- 182,500 equipment logbook entries per year,
- 365,000 equipment logbook entries after 2 years, and
- an average of 73 update-related equipment logbook entries per entity after 2 years.

Performance of Unit Class Export 1 (SR1095.10.2.1)

Based on Reference scenario 1 (page 6), the export of 1 unit class does not take longer than 5 seconds.

Performance of Unit Class Import 1 (SR1095.10.2.2)

Based on Reference scenario 1 (page 6), the import of 1 unit class does not take longer than 5 seconds for an initial import.

Performance of Filter Class Export 1 (SR1095.10.2.5)

Based on Reference scenario 1 (page 6), the export of 1 filter class does not take longer than 15 seconds.

Performance of Filter Class Import 1 (SR1095.10.2.6)

Based on Reference scenario 1 (page 6), the import of 1 filter class does not take longer than 15 seconds for an initial import.

Recipe Execution - Performance (SR1095.10.7+)

For details about the execution with PharmaSuite, see "Functional Requirement Specification Execution Framework" [A3] (page 39).

Reference scenario 1

The reference recipe holds a total of:

- 1 Dispense unit procedure,
- 1 Dispense operation
- All Dispense-related phases (D Identify material, D Select scale, D Tare, D Weigh, D Release scale, D Print report)
- 10 materials to be dispensed

Boundary scenario 1

The boundary recipe holds a total of:

- up to 30 unit procedures,
- up to 35 operations,
- up to 3,000 phases, and
- up to 12,000 process parameters.

Its procedural structure includes various branches and loops.

Performance of Recipe Execution 1 (SR1095.10.7.1)

Based on the Boundary scenario 1 (page 8) and phase completion through automated test, the related control recipe is executable within 60 minutes.

Printing of Barcode Label and Dispensing Report - Performance

For details about the barcode labels, see "Functional Requirement Specification Runtime Data Management" [A8] (page 39), "Functional Requirement Specification Dispense and Inline Weighing" [A9] (page 39), and "Functional Requirement Specification Output Weighing" [A10] (page 39).

For details about the dispensing report, see "Functional Requirement Specification Dispense and Inline Weighing" [A9] (page 39).

Performance of Printing (SR1090.4.6)

Printing a barcode label does not take longer than 3 seconds (on average) and less than 5 sec (in 80% of the cases).

A reference Dispensing report contains of four pages and the reference printer prints at least 15 pages per minute.

Based on the reference report, printing the Dispensing report does not take longer than 30 seconds (in the average) and less than 60 seconds (in 80% of the cases).

EBR Server - Performance (SR1095.10.3+)

For details about the execution with PharmaSuite, see "Functional Requirement Specification Execution Framework" [A3] (page 39).

Reference scenario 1

For the performance test of the EBR server, the reference scenario is as follows:

- Size of recipe:
 - 4 unit procedures (2 parallel)
 - Each with 4 operations (2 parallel)
- Order explosion and order execution take place in parallel
- 50 Production Execution Clients
- 500 orders (10 sequential orders per client)

The focus of the performance test is on the behavior of the EBR server. Therefore, the execution of operations is performed as a virtual execution on a virtual client, which results in a zero second response time from a client.

Performance of Order Processing 1 (SR1095.10.3.1)

Based on Reference scenario 1 (page 9), the EBR server is capable to process 500 orders in less than 50 minutes (10 orders or 160 operations per minute).

TOM Server - Performance (SR1095.10.5+)

For details about the execution with PharmaSuite, see "Functional Requirement Specification Execution Framework" [A3] (page 39).

Reference scenario 1

For the performance test of the TOM server, the reference scenario is as follows:

- Test case includes:
 - 3 times start and removal of ETO templates
 - 10 times start and completion of ETO runs
 - 10 times escalation and expiration of ETO runs
- Scale-up consists of:
 - 50 Production Execution Clients running in parallel, all connected to the same TOM server
 - 10 consecutive runs of the test case on each client

The reference scenario results in a total of

- 1,500 ETO templates,
- 5,000 completed ETO runs, and
- 5,000 expired ETO runs,

that are managed by the TOM server.

The focus of the performance test is on the behavior of the TOM server. Therefore, the execution of operations is performed as a virtual execution on a virtual client, which results in a zero second response time from a client.

Performance of ETO Processing 1 (SR1095.10.5.1)

Based on Reference scenario 1 (page 10), the TOM server is capable to process ETO templates and runs in less than 30 minutes (50 templates per minute, 330 completed and expired runs per minute).

Exception Dashboard - Performance (SR1095.10.8+)

For details about the Exception Dashboard, see "Functional Requirement Specification Review and Approval" [A4] (page 39).

Reference scenario 1

The reference recipe holds a total of:

- 3 unit procedures,
- 5 operations,
- 20 phases, and
- 100 process parameters.

Performance of Review by Exception 1 (SR1095.10.8.1)

Based on the Reference scenario 1 (page 11) and a finished order, opening a unit procedure with 50 exceptions within the Details panel does not take longer than 50 seconds.

Performance of Batch Report Printing 1 (SR1095.10.8.2)

Based on the Reference scenario 1 (page 11), the response time for printing the Batch Report of a finished order (time between "trigger the print function" and "display of the preview") does not take longer than 1 minute.

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Compliance Requirements (SR1095.50+)

The detailed list of the PharmaSuite compliance requirements is structured in the following sections:

- **Audit Trail (SR1095.50.1)** requirements (page [13](#))
- **Electronic Signatures (SR1095.50.2)** requirements (page [13](#))
- **Access Control (SR1095.50.3)** requirements (page [16](#))

Audit Trail (SR1095.50.1)

The system shall provide functionality for audit trail.

The following data shall be included as part of any human-readable form of an audit trail record:

- Transaction type (insert, update, delete).
- Timestamp.
- Logged-in user (login name, first and last name).
- Link to electronic signature(s), if applicable
- Transaction type-specific data (shall be defined for each transaction type)

Audit trail data is available for:

- Users, **Audit Trail (SR1095.50.3.2.1)** requirement (page [18](#))
- User groups, **Audit Trail (SR1095.50.3.3.1)** requirement (page [20](#))
- Access privileges, **Audit Trail (SR1095.50.3.4.1)** requirement (page [20](#))
- Materials, **Audit Trail (SR3146.1.1.5)** requirement in "Functional Requirement Specification Data Management" [A2] (page [39](#))

Electronic Signatures (SR1095.50.2)

The system shall support electronic signatures.

For electronics signature in the context of EBR and the Production Execution Client, see **Signatures for Exceptions and Comments (SR1079.2.1)** and **Signature for Phase Completion (SR1089.4+)** in "Functional Requirement Specification Execution Framework" [A3] (page [39](#)).

Single and Double Electronic Signatures (SR1095.50.2.1)

The system shall support single and double electronic signatures. Double electronic signatures support the witness role.

When an electronic signature is required, the system shall request the user to authenticate with his user name and password. In case of a double electronic signature, two consecutive electronic signatures are required from two different users.

Whether a single or a double electronic signature is required shall be defined in the system configuration or during recipe or workflow definition.

The screenshot shows a single 'Electronic Signature' dialog box. It has a title bar with a menu icon and a close button. The dialog contains five labeled text input fields: 'Description' (containing 'Remove object locks'), 'Reason' (containing 'Confirmation'), 'User Name' (containing 'dgerr'), 'Password' (a yellow masked field), and 'Comment' (empty). At the bottom are 'OK' and 'Cancel' buttons. The identifier 'Identifier:PMC_ES_Remove...' is displayed at the bottom right.

Figure 1: Example of a single electronic signature

The screenshot shows two 'Electronic Signature' dialog boxes side-by-side, representing a double signature process. Both have the same title bar. The left dialog has 'Description: Change batch attributes', 'Reason: Confirmation', 'User Name: dgerr', a yellow 'Password' field, and an empty 'Comment' field. The right dialog has 'Description: Change batch attributes', 'Reason: Approval', a yellow 'User Name' field, a yellow 'Password' field, and an empty 'Comment' field. Both have 'OK' and 'Cancel' buttons. The identifier 'Identifier:Sig_ChangeBatchAttributes' is shown at the bottom right of the right dialog.

Figure 2: Example of a double electronic signature to support witness role

For electronic signatures, the authentication requires the login name and password of a user as used for system login.

COMMENTS (SR1095.50.2.1.1)

The system shall allow to enter a comment along with an electronic signature.
This does not apply to phase completion signatures during execution in the Production Execution Client.

Whether a comment is mandatory, optional, or not supported shall be defined in the system configuration or during recipe or workflow definition.

LINKAGE TO ELECTRONIC RECORD (SR1095.50.2.1.2)

The system shall ensure that the linkage of an executed electronic signature to the electronic records for which the signature was given is unambiguous and incorruptible.

Configuring Electronic Signatures (SR1095.50.2.2)

The system shall allow the system administrator and system integrator to configure the following parameters of an electronic signature:

- Signature type (single, double)
- Meaning of the signature (description)
- Required user role of the user who signs (reason)
- Comment (no comment, optional, mandatory)
This does not apply to phase building blocks which are executed in the Production Execution Client.
- Entitled user group
- Signature level (disabled, implicit w/o warning, explicit, warning w/ implicit, error (= no way to proceed))
This does not apply to phase building blocks which are executed in the Production Execution Client.
- With the **implicit w/o warning** level, the system supports silent electronic signatures.
If configured, instead of requesting an electronic signature, the system shall record the logged-in user automatically.

Performing Electronic Signatures (SR1095.50.2.3)

When the system requests an electronic signature it shall unambiguously provide the signer with the following information:

- Meaning of the signature
- Required user role (e.g. Confirmation, Approval)
- Scope of the signature

- Comment (as defined for the signature):
 1. no comment
 2. optional comment
 3. mandatory comment
- In case the user is the second signer of a double signature, the data entered by the first signer must be visible for the second signer.

Data of an Electronic Signature (SR1095.50.2.4)

The system shall add the following data of a given electronic signature (single or double) to any human-readable form of the affected electronic record (e.g. printout, electronic display):

- Meaning of the signature
- Role of the user who signed
- Scope of the signature as seen by the signer (context data related to the signature)
- Signer, timestamp, and comment (if mandatory or provided by the user)

Access Control (SR1095.50.3)

The system shall provide functionality for access control.

For access control requirements of specific system parts (e.g. access rights for editing or viewing data in Data Manager), please refer to

- "Functional Requirement Specification Execution Framework" [A3] (page 39)
- "Functional Requirement Specification Runtime Data Management" [A8] (page 39)
- "Functional Requirement Specification Review and Approval" [A4] (page 39)
- "Functional Requirement Specification Recipe and Workflow Management" [A1] (page 39)
- "Functional Requirement Specification Data Management" [A2] (page 39)

Authentication for System Login and Electronic Signature (SR1095.50.3.1)

The system shall ensure that only users with sufficient privileges have access to the system and are allowed to perform an electronic signature.

VERIFICATION DURING AUTHENTICATION (SR1095.50.3.1.1)

The system shall perform the following verification checks:

- Login name belongs to a valid user account.

- Password is correct.
- Password has not expired.
- User account is activated.
- Maximum number of invalid authentication attempts is not exceeded.
- User has sufficient access privileges to access the application or perform the electronic signature.

INVALID AUTHENTICATION ATTEMPT (SR1095.50.3.1.2)

In case of an invalid authentication attempt, the system shall display system messages to the user that do not give any hints to a possible intruder.

LOGGING UNSUCCESSFUL AUTHENTICATION ATTEMPTS (SR1095.50.3.1.3)

The system shall log unsuccessful authentication (login and electronic signature) attempts.

It shall be possible to notify a dedicated responsible person (or persons) of the event. The configurable notification shall include the following information:

- Used login name (identification code).
- Change of user account status (when the number of invalid login attempts is reached).
- Description of unsuccessful authentication.
Example: Attempt to authenticate with a locked account.

The used medium for notifications shall be e-mail.

NUMBER OF INVALID AUTHENTICATION ATTEMPTS EXCEEDED (SR1095.50.3.1.4)

In case the number of invalid authentication attempts of a user account is exceeded, any further authentication request shall be blocked by the system. The user must contact the system administrator to reset the number of invalid authentication attempts.

MULTIPLE LOGINS PER USER (SR1095.50.3.1.5)

The system shall allow that one user is logged in in more than one session at the same moment.

VIEWING SYSTEM ACCESS HISTORY (SR1095.50.3.1.6)

The system shall track system access of users.

It allows a user to view the users with the login (successful, invalid) and logout information. The information includes (but is not limited to) the following properties:

- Timestamp
- Login name

- Transaction type (login, logout, failed login)
- Database name

Management of Users (SR1095.50.3.2)

The system allows the system administrator to define and manage users and their access privileges.

AUDIT TRAIL (SR1095.50.3.2.1)

The system allows a user to view audit trail data recorded for users.

DEFAULT USER ACCOUNT SETTINGS (SR1095.50.3.2.2)

The system shall allow the system administrator to define default settings for user accounts:

- Maximum number of invalid authentication attempts: [0..n]
Default: 5
- Validity timeframe for user accounts: [0..n]
Default: 365 days
- Expiration time for passwords: [0 days..n years]
Default: 365 days
- Notification of password expiration at login: [0 days..n days]
Default: 10 days
- Password strength:
 - Minimum length
Default: 6
 - Maximum length
Default: 20
 - Minimum number of lowercase characters
Minimum number of uppercase characters
Minimum number of numeric characters
Minimum number of special characters
Default: 1..n [a..z] + 0..n [A..Z] + 1..n [numeric] + 1..n [special character]
 - Minimum number of previous password before reuse: [0..n]
Default: 3

UNIQUE LOGIN NAME (SR1095.50.3.2.3)

The system shall enforce the system administrator to define a unique login name when creating a user account.

With a unique login name (identification code), the combination of it and a password is always unique.

IMMUTABLE LOGIN NAME (SR1095.50.3.2.4)

The system shall prevent that a login name of a user account can be changed after the account has been created.

VALIDITY TIMEFRAME (USER ACCOUNT) (SR1095.50.3.2.5)

The system shall allow the system administrator to define the date when a user account becomes invalid.

The system shall use default values for the end date of the validity timeframe as defined in the system configuration according to the **Default User Account Settings (SR1095.50.3.2.2)** requirement (page 18).

ACTIVITY STATUS (USER ACCOUNT) (SR1095.50.3.2.6)

The system shall allow the system administrator to activate and deactivate a user account.

NUMBER OF INVALID AUTHENTICATION ATTEMPTS (SR1095.50.3.2.7)

The system shall allow the system administrator to define how many false authentication attempts are allowed for a user account before it is locked.

The system shall use default values for the maximum number of invalid authentication attempts as defined in the system configuration according to the **Default User Account Settings (SR1095.50.3.2.2)** requirement (page 18).

PASSWORD STRENGTH (SR1095.50.3.2.8)

The system shall allow the system administrator to define the required password strength. The system shall use the default values as defined in the system configuration according to the **Default User Account Settings (SR1095.50.3.2.2)** requirement (page 18).

EXPIRY DATE (PASSWORD) (SR1095.50.3.2.9)

The system shall allow the system administrator to define an expiry date for the password of a user account.

The system shall use the rules for generating the password expiry date as defined in the system configuration according to the **Default User Account Settings (SR1095.50.3.2.2)** requirement (page 18).

ASSIGNING USERS TO USER GROUPS (SR1095.50.3.2.10)

The system shall allow the system administrator to assign a user to one or multiple user groups and to revoke the assignment.

The system shall not limit the number of user groups to which a user can be assigned. A user shall inherit the access privileges of the user group to which the user is assigned.

Management of User Groups (SR1095.50.3.3)

The system allows the system administrator to define and manage user groups and their access privileges.

AUDIT TRAIL (SR1095.50.3.3.1)

The system allows a user to view audit trail data recorded for user groups.

Management of Access Privileges (SR1095.50.3.4)

The system allows the system administrator to define and manage access privileges (access rights and electronic signatures).

AUDIT TRAIL (SR1095.50.3.4.1)

The system allows a user to view audit trail data recorded for access privileges.

ACCESS PRIVILEGES FOR PRE-DEFINED SET OF FUNCTIONS (SR1095.50.3.4.2)

The system shall provide access privileges for a set of pre-defined functions.
This set of functions shall be extendable by a system integrator.

PRE-DEFINED ROLES AND PROFILES (SR1095.50.3.4.3)

The system shall provide a set of pre-defined roles, mapped to user groups, with pre-defined access privileges according their assigned profiles.
This set of user groups and access privileges shall be changeable and extendable by a system integrator.

VIEWING ACCESS PRIVILEGES OF USERS OR USER GROUPS (SR1095.50.3.4.4)

The system allows a user to view the access privileges of a user or user group at any defined moment.

This moment can be in the past until the present time.

VIEWING USERS OR USER GROUPS OF ACCESS PRIVILEGES (SR1095.50.3.4.5)

The system allows a user to view the users or user groups of an access privilege at any defined moment.

This moment can be in the past until the present time.

Password Security (SR1095.50.3.5)

The system shall ensure confidentiality of passwords that are managed in the system.
Only the related person shall know his personal password.

An exception is the initial password.

PASSWORD ENCRYPTION (SR1095.50.3.5.1)

The system shall fulfill the following items related to password encryption:

- Passwords are always stored encrypted.
- The system must not provide a function to decrypt an encrypted password. The encryption is a one-way encryption.
- The system encrypts passwords that need to be sent over a network that can be accessed by other users.

PASSWORD MUST NOT BE DISPLAYED (SR1095.50.3.5.2)

The system shall never display a password in plain text on any kind of user interface including printouts.

Password Change (SR1095.50.3.6)

The system shall support that the password can be changed.

PASSWORD CHANGE - USER-INITIATED (SR1095.50.3.6.1)

The system shall provide functionality for the user to change his password. The system requests the user to enter his old password and the new password twice. The new password shall fulfill the password requirements as defined in the system configuration according to the **Default User Account Settings (SR1095.50.3.2.2)** requirement (page 18).

PASSWORD CHANGE - AFTER LOGIN (SR1095.50.3.6.2)

The system shall allow a user to change the password directly after a successful login. This is important in case the password is going to expire within the next n days and the user needs to have the system available immediately.

PASSWORD CHANGE - SYSTEM-FORCED (SR1095.50.3.6.3)

In case the password has expired, the system shall force the user to change his password at login. There shall be no possibility to access the system without changing the expired password.

The system requests the user to enter his old password and the new password twice. The new password shall fulfill the password requirements as defined in the system configuration according to the **Default User Account Settings (SR1095.50.3.2.2)** requirement (page 18).

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Internationalization and Localization (SR1095.1+)

PharmaSuite is marketed globally. This requires the software to support internationalization (I18N) and localization (L10N).

Internationalization (I18N) is the process of planning and implementing products and services so that they can easily be adapted to specific local languages and cultures, a process called localization. The internationalization process is sometimes called translation or localization enablement.

Localization (L10N) is a process of adapting a product or service to a particular language, culture, and desired local "look-and-feel".

The requirements specific to internationalization and localization are structured in the following sections:

- **Skills for Localization (SR1095.1.1)** requirements (page [23](#))
- **Locale-related** requirements (page [24](#))

Skills for Localization (SR1095.1.1)

The system shall be localizable into a specific language or terminology by a system integrator.

It shall be possible that a person used to work with office software without knowledge of the system or specific IT knowledge is able to modify choice lists contents and text labels of the user interface and reports. It shall be possible for a system integrator to export this information from the system. The result of the modification must be in a format that it can be integrated to the system by a system integrator.

Tool Support (SR1095.2.4)

The system shall provide tools that allow a system integrator to localize the user interface, reports, barcode labels, system messages, etc. for a particular language or terminology.

Internationalized Output Media (SR1095.2)

The system shall provide internationalization for those parts of the system that are presented to common users of the system, so that it is possible to localize them into a particular language or terminology. The user roles do not include system integrators, project engineers, etc.

The internationalized parts include the user interface, online help, reports, barcode labels, system messages, installation software, and user manuals. Internationalization is not limited to text but also includes data formats like date, time, and number formats.

Locale-related Requirements

A locale defines the regional settings with specific parameters. For example language of the user interface, keyboard layout, number format, date and time format, paper size.

Supported Locales (SR1095.1.2)

The system shall support localization into the following locales (sequence not related to priority):

- English (United States) - Default
- German (Germany)
- Spanish (Spain-Modern Sort)
- Italian (Italy)
- French (France)
- Chinese (PRC)
- Korean
- Japanese
- Russian (Russia)
- Hebrew
- Hindi

Provided Locales (SR1095.2.3+)

English (United States) (SR1095.2.3.1)

The default locale of the system shall be English (United States):

- Default date format: mm/dd/yyyy
- Default digit grouping symbol: , (comma)
- Default decimal separator: . (dot)
- Default time format: hh:mm:ss (AM/PM)

Locales within One Installation (SR1095.2.1)

The system shall support multiple locales within one installation.

Choice of Locale (SR1095.2.2)

The system shall allow the user to choose the locale to be used at when logging on to the system.

Configurable Locale-specific Settings (SR1095.1.3)

The system shall allow the configuration of the following locale-specific settings:

- Date format
- Number format (digit grouping symbol and decimal separator)
- Time format

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Installation, Backup, and Upgrade (SR1095.4+)

➤ For recent changes, see revision history (page 41).

The PharmaSuite Express Edition Installer significantly simplifies the installation process for smaller systems. The goal is to enable a system administrator or integrator to install the software without the need of manual configuration and rework by guiding him with the installation program through the required installation steps. An initial installation of a system should not take longer than 1 hour.

The following limitations apply:

- Up to the current release, the PharmaSuite Express Edition Installer is limited to the initial installation of PharmaSuite. The installation of maintenance releases, minor software updates, or major software upgrades into an existing PharmaSuite system will be covered by the upgrade installer, which has been introduced with PharmaSuite 4.0.
- There will be no PharmaSuite Enterprise Edition Installer. Instead, the installation will be into an existing FactoryTalk ProductionCentre system according to the "Technical Manual Installation - Enterprise Edition" of PharmaSuite.

Installation of PharmaSuite (SR1095.4.1)

The system shall allow to be installed by system administrators.

Initial Installation (SR1095.4.1.1)

The installation of the system shall support initial installations (FactoryTalk ProductionCentre platform or applications not yet installed).

Existing Database (SR1095.4.1.7)

The installation shall support installation into an existing database installation (other instances already installed on the database).

Installation Result (SR1095.4.1.8)

The installation shall inform the user if the installation was successful or not.

LOGGING OF INSTALLATION (SR1095.4.1.8.1)

The installation shall write all installation steps into a log file including the status of each step and, if applicable, the reason for failure.

Scope of Express Edition Installer (SR1095.4.1.6)

For Express Edition, an installation program shall cover all components of the system necessary to start operation. This includes:

- FactoryTalk ProductionCentre platform (through included FactoryTalk ProductionCentre installer)
- Application server (through included FactoryTalk ProductionCentre installer)
- PharmaSuite applications and capabilities (PharmaSuite software and configuration data)

CONFIGURATION PARAMETERS OF EXPRESS EDITION INSTALLER (SR1095.4.1.9)

For Express Edition, the installation shall prompt the administrator for all parameters required for the installation, like database server name and password, database name, etc.

Installation Manual (SR1095.4.1.4)

The system shall provide an installation manual that supports a system administrator when installing the system.

Backup (SR1075.6)

- For recent changes, see revision history (page 41).

The system shall allow to backup data and system.

Backup of a Running System (SR1075.6.1)

- For recent changes, see revision history (page 41).

The system shall allow to backup the data stored in the database during production time with an appropriate third-party backup tool. The backup does not include data in the local client's cache.

Backup and Restoring (SR1075.6.3)

- For recent changes, see revision history (page 41).

It shall be possible to completely backup and restore both, system and data, with a third-party tool. The configuration and execution of the backup and restore environment is done by IT subject matter experts.

Upgrade Tool (SR1075.8+)

- For recent changes, see revision history (page 41).

The system shall provide an upgrade tool to update the system and to migrate the data.

Basic Upgrade Engine (SR1075.8.1)

- For recent changes, see revision history (page 41).

The basic upgrade engine allows the upgrade of a PharmaSuite releases in five stages:

- Stage 1 verifies the completeness of the upgrade package.
- Stage 2 determines the planned changes required to upgrade from the base to the new version.
- Stage 3 checks the planned changes against the target system and reports non-resolvable conflicts as well as potential conflicts (based on project extensions).
- Stage 4 performs the system update and the data migration on the target system based on the results of stage 3.
- Stage 5 verifies the results of the upgrade process.

Upgrade - Continue execution (SR1075.8.3.5)

- For recent changes, see revision history (page 41).

Along with the upgrade of the system to a following new software version, the system shall allow to execute and finish orders with the new release that originally have been started with the previous release.

System Update Tasks (SR1075.8.2+)

- For recent changes, see revision history (page 41).

The upgrade tool shall provide the following system update tasks:

Update task for DSX objects (generic DSX) (SR1075.8.2.1)

- For recent changes, see revision history (page 41).

The upgrade tool shall support an update task for objects that are part of FactoryTalk ProductionCentre DSX format (generic DSX).

Update task for Message objects (SR1075.8.2.2)

- For recent changes, see revision history (page 41).

The upgrade tool shall support an update task for message packs. The conflict handling shall be done on message ID level.

Update task for EAR files (SR1075.8.2.3)

- For recent changes, see revision history (page 41).

The upgrade tool shall support an update task for EAR artifacts.

Update task for ATRow objects (SR1075.8.2.4)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for ATRow artifacts.

Update task for Library objects (SR1075.8.2.5)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for library artifacts.

Update task for Image objects (SR1075.8.2.6)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for image artifacts.

Update task for FSM (Flexible State Model) objects (SR1075.8.2.7)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for FSM artifacts.

Update task for AT Definition objects (SR1075.8.2.8)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for AT Definition artifacts.

Update task for UDA Definition objects (SR1075.8.2.9)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for UDA Definition artifacts.

Update task to report customer notifications (SR1075.8.2.11)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task to report customer modifications.
This task shall only be applicable to **Stage 3** (Checking Target System for Conflicts) and create a report with the list of objects that have been added to, modified in, or deleted from the customer system.

Update task for Application objects (SR1075.8.2.12)

- For recent changes, see revision history (page 41).
- The upgrade tool shall support an update task for Application artifacts.

Update task for Versioning Configuration objects (SR1075.8.2.13)

- For recent changes, see revision history (page 41).

The upgrade tool shall support an update task for Versioning Configuration artifacts.

Update task for Report Design objects (SR1075.8.2.14)

- For recent changes, see revision history (page 41).

The upgrade tool shall support an update task for Report Design artifacts.

Update task for Access Privilege objects (SR1075.8.2.15)

- For recent changes, see revision history (page 41).

The upgrade tool shall support an update task for Access Privilege artifacts.

System Update and Data Migration Packages (SR1075.8.3+)

- For recent changes, see revision history (page 41).

The upgrade tool shall provide system update and data migration packages specific to a PharmaSuite version.

System update from the previous software version (SR1075.8.3.25)

- For recent changes, see revision history (page 41).

The system update package shall allow to update a PharmaSuite system from the previous software version to the subsequent software version by applying functional and structural changes and adding new functions.

For an overview of the default system update tasks, please refer to the **System Update Tasks (SR1075.8.2+)** requirement (page 29).

For a system update from PharmaSuite 8.3 to PharmaSuite 8.4, the following specific system update task is available:

- Migrate the **S88EquipmentClassStatusGraph** flexible state model

Data migration from the previous software version (SR1075.8.3.26)

- For recent changes, see revision history (page 41).

The data migration package shall allow to migrate specific data that already resided in a previous software version of PharmaSuite to be suitable for use with the new or changed system functions of the subsequent software version.

For a data migration from PharmaSuite 8.3 to PharmaSuite 8.4, the following specific data migration tasks are available:

- Migrate transition conditions of activity sets
- Update usage list-related tables

Extensibility (SR1075.8.5+)

- For recent changes, see revision history (page [41](#)).

Upgrade packages (SR1075.8.5.1)

- For recent changes, see revision history (page [41](#)).

The system upgrade tool shall allow a system integrator to adapt system update and data migration packages, so that customized releases of PharmaSuite can be upgraded with the extended upgrade packages.

This extension use case shall be documented and informally tested.

Pre-install Configuration (SR3146.7)

The system shall support specific pre-install configurations by a system integrator.

User Experience (SR1095.5+)

Messages (SR1095.5.1.1)

The system shall display only meaningful system and error messages to the user. Pure technical messages shall be avoided.

Example: The operator is generally not interested that a unique constraint in a database table is violated, but that the batch identifier is not unique.

Device-specific Usability (SR1095.5.1.4+)

Touch Screen (SR1095.5.1.4.2)

The system shall provide a user interface for the Production Execution Client that allows to operate the system via touch screen.

Virtual Keyboard (SR1095.5.1.4.2.1)

As an EBR operator I need to be able to enter data via a virtual keyboard ("soft keyboard"), so that I can operate the system via touch screen only.

Look and Feel (SR1095.5.7)

The Look and Feel of the system shall be state-of-the-art.

GUI Design (SR1095.5.7.3)

The GUI shall be designed in a way that supports easy readability of key information, in support of the special needs of the operator. This includes structure, design and sizing of forms and panels.

GUI Adaptation (SR1095.5.7.1)

The user interface of PharmaSuite shall be adaptable by simple means, allowing to adjust colors, fonts, and font sizes.

Dialog Boxes (SR1095.5.7.2)

The dialog boxes used within PharmaSuite shall all have a common appearance in terms of colors, fonts, and button style.

Grid Components (SR1095.5.7.4)

The system shall support state-of-the-art grid components. When functionally suitable, this includes:

- sorting by rows,
- manual adjustment of sequence of columns, and
- support of multi-selection (per configuration).

Extensibility (SR1095.5.1.2+)

System - Customization (SR1095.5.1.2.1)

The system shall be customizable in a way that a system integrator is able to:

- add or modify workflows,
- add, remove, or modify data object attributes,
- add, remove, or modify fields and labels of the user interface,
- rename labels of the user interface, and
- add, remove, or modify reports and barcode labels.

It must be possible to perform the modifications "offline" in a development system and to integrate them into an existing customer system. The respective guidelines shall be documented in Technical Manuals.

System - View Model-related Objects (SR1095.5.1.2.3)

The system shall allow a system integrator to add attributes to a view model-related object without the need of Java coding.

This extension use case shall be documented and tested.

System - Material Attribute (SR1095.5.1.2.9)

The system shall allow a system integrator to add a material/batch attribute across the system and to implement business logic that makes use of this attribute.

This extension use case shall be documented and informally tested. Testing shall be based on the example of a **Second potency** attribute.

System - ERP BOM Attribute (SR1095.5.1.2.11)

The system shall allow a system integrator to add a new attribute to the ERP BOM, which is automatically populated into the related new attribute of the material parameter.

In case the same attribute also exists as material master data, the default value for the material parameter will be taken from the material master data in case no value exists in the ERP BOM.

This extension use case shall be documented and informally tested. Testing shall be based on the example of a **Second potency** attribute.

System - Batch Status (SR1095.5.1.2.10)

The system shall allow a system integrator to add a batch status to the flexible state model (FSM) and to implement business logic that makes use of the new status.

This extension use case shall be documented and informally tested. Testing shall be based on the example of a **Conditionally released** status of a batch.

Master Recipe - Approval Workflow (SR1095.5.1.2.13)

The system shall allow a system integrator to adapt the approval workflow of a master recipe according to project needs. This includes the ability to

- add a new status that has the same semantic meaning as one of the given **Edit** and **Verification** statuses,
- adjust status transitions accordingly, and
- rename statuses and status transitions (actions).

This extension use case shall be documented and informally tested.

Expression Editor - Context Data (SR1095.5.1.2.14)

The system shall allow a system integrator to enhance the context data-related capability of the Expression editor. This includes the ability to

- add a new attribute to the list of available attributes for a given object-related function and
- add a new function to access the context data of a new object.

This extension use case shall be documented and informally tested.

Building Block - Name Spaces (SR1095.5.1.2.4)

The system shall support the deployment of building blocks from different system integrators by definition of name spaces.

This extension use case shall be documented and tested.

Building Block - Mainstream GUI Programming (SR1095.5.1.2.6)

For the implementation of a GUI of a System Building Block, the system shall support mainstream Java GUI programming mechanisms, rather than being dependent on proprietary components.

Execution - Autostart (SR1095.5.1.2.2)

For a system integrator, the system shall allow to configure specific activity sets that are started automatically upon start of the Production Execution Client. One activity set can

be started multiple times in case multiple instances are allowed for this activity set (configurable).

Execution - Select Order Grid (SR1095.5.1.2.5)

For the **Select order** grid of the Batch, Device, or Workflow Processing actions of the Cockpit, the system shall allow a system integrator to add project-specific fields/columns to the grid and to adjust the pre-defined filter criteria, according to project-specific business rules.

This extension use case shall be documented and tested.

Execution - EBR Header Bar (SR1095.5.1.2.8)

The system shall allow a system integrator to change the content of the EBR header bar of the Production Execution Client.

Reporting (SR1095.5.1.2.16+)

MASTER RECIPE REPORT VARIANTS (SR1095.5.1.2.16.1)

The system shall allow a system integrator to create and print a master recipe report variant from Recipe Designer in addition to the standard master recipe report.

This extension use case shall be documented and informally tested.

MASTER WORKFLOW REPORT VARIANTS (SR1095.5.1.2.16.2)

The system shall allow a system integrator to create and print a master workflow report variant from Workflow Designer in addition to the standard master workflow report.

This extension use case shall be documented and informally tested.

BATCH REPORT VARIANTS (SR1095.5.1.2.16.3)

The system shall allow a system integrator to create and print a batch report variant from the Production Management Client and the Production Response Client in addition to the standard batch report.

This extension use case shall be documented and informally tested.

WORKFLOW REPORT VARIANTS (SR1095.5.1.2.16.4)

The system shall allow a system integrator to create and print a workflow report variant from the Production Management Client in addition to the standard workflow report.

This extension use case shall be documented and informally tested.

Precision of Numerical Data (SR1095.5.1.3+)

Fixed Point Numbers (SR1095.5.1.3.1)

The system shall support decimal fixed-point numbers for the definition of input and output parameters to support the control of scaling of decimal values.

Floating Point Arithmetic (SR1095.5.1.3.2)

For calculations, the system shall use decimal floating point arithmetic that is compliant to the IEEE 754R standard.

Errors at Calculation (SR1095.5.1.3.3)

Errors at calculation (division by zero, overflow, underflow, inexact, invalid operation), as defined in the IEEE 754R standard, shall be reported by the system.

Input Data Errors (SR1095.5.1.3.4)

The system shall not consider input data errors at calculation, i.e. the scale of a calculation result is not adapted by the system. The correct setting of the scale of a calculation result lies within the responsibility of the system user.

Rounding (SR1095.5.1.3.5)

At rounding, the following method shall apply (round-half-up): if the discarded digits are greater than or equal to a half (0.5) of a one in the next left position then the result coefficient shall be incremented by 1 (rounded up) - otherwise the discarded digits are ignored.

Concurrent Data Access - Pessimistic Locking (SR1095.5.1.5)

Concurrent data access shall not lead to data inconsistency.
The system shall support pessimistic locking mechanisms.

Example:

Record 1 is locked by user 1; user 1 has read and write rights.

For user 2, who comes later than user 1, the system assures that user 2 cannot modify or delete the locked record.

In the Production Management Client and Recipe and Workflow Designer, the system still allows user 2 to read the record while it is locked by user 1.

Restart and Recovery (SR1095.8)

The system shall support "save points".

At a save point, the system generates consistent data in the system accessible by other users and functions.

The default save points are defined in the application. They are not configured nor shall be configurable. The locations of the save points are defined in the requirements of the respective pre-defined workflows. In case of a system breakdown, all modifications done up to there are rolled back if no save point has been reached in the respective workflow.

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Reference Documents

The following documents are available from the Rockwell Automation Download Site.

No.	Document Title	Part Number
A1	PharmaSuite Functional Requirement Specification Recipe and Workflow Management	PSFRSRD-RM008E-EN-E
A2	PharmaSuite Functional Requirement Specification Data Management	PSFRSDM-RM004E-EN-E
A3	PharmaSuite Functional Requirement Specification Execution Framework	PSFRSEF-RM004E-EN-E
A4	PharmaSuite Functional Requirement Specification Review and Approval	PSFRSRA-RM004E-EN-E
A5	PharmaSuite Functional Requirement Specification Equipment Automation Phases	PSFRSEA-RM004E-EN-E
A6	PharmaSuite Technical Manual Installation - Enterprise Edition	PSEN-IN008E-EN-E
A7	PharmaSuite Technical Manual Installation - Upgrade	PSUP-IN005E-EN-E
A8	PharmaSuite Functional Requirement Specification Runtime Data Management	PSFRSRT-RM004E-EN-E
A9	PharmaSuite Functional Requirement Specification Dispense and Inline Weighing	PSFRSDI-RM006E-EN-E
A10	PharmaSuite Functional Requirement Specification Output Weighing	PSFRSOW-RM002E-EN-E

TIP

To access the Rockwell Automation Download Site, you need to acquire a user account from Rockwell Automation Sales or Support.

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

Document Information

The document information covers various data related to the document.

Approval

This document has been approved electronically via the Rockwell Automation Document Management System (DMS). The required approvers of this document include the following:

Name	Role
Martin Dittmer	Product Manager
Steffen Landes	Development Manager
Martin Irmisch	Test Manager

In addition, the electronic document approval via DMS is confirmed by a handwritten signature of all approvers in the Quality Document when the release is completed. The Quality Document summarizes the quality-related planning activities and results of a PharmaSuite release.

Version Information

Object	Version
PharmaSuite	8.4
Functional Requirement Specification	1.0

Revision History

The following table describes the history of this document.

Changes related to the document:

Object	Description	Document
Backup (SR1075.6) (page 28)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Backup of a Running System (SR1075.6.1) (page 28)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Backup and Restoring (SR1075.6.3) (page 28)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Upgrade Tool (SR1075.8+) (page 28)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Basic Upgrade Engine (SR1075.8.1) (page 29)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
System Update Tasks (SR1075.8.2+) (page 29)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for DSX Objects (Generic DSX) (SR1075.8.2.1) (page 29)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Message Objects (SR1075.8.2.2) (page 29)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for EAR Files (SR1075.8.2.3) (page 29)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for ATRow Objects (SR1075.8.2.4) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Library Objects (SR1075.8.2.5) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Image Objects (SR1075.8.2.6) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for FSM (Flexible State Model) Objects (SR1075.8.2.7) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0

Object	Description	Document
Update Task for AT Definition Objects (SR1075.8.2.8) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for UDA Definition Objects (SR1075.8.2.9) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task to Report Customer Notifications (SR1075.8.2.11) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Application Objects (SR1075.8.2.12) (page 30)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Versioning Configuration Objects (SR1075.8.2.13) (page 31)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Report Design Objects (SR1075.8.2.14) (page 31)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Update Task for Access Privilege Objects (SR1075.8.2.15) (page 31)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
System Update and Data Migration Packages (SR1075.8.3+) (page 31)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management". Rephrased to be more precise about system update and data migration.	1.0
Upgrade - Continue Execution (SR1075.8.3.5) (page 29)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Extensibility (SR1075.8.5+) (page 32)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0
Upgrade Packages (SR1075.8.5.1) (page 32)	Editorial Moved from "PharmaSuite 8.0 - DRS Information Management".	1.0

Changes related to "Performance Requirements" (page 3):

Object	Description	Document
--	---	---

Changes related to "Compliance Requirements" (page 13):

Object	Description	Document
---	---	---

Changes related to "Internationalization and Localization" (page 23):

Object	Description	Document
---	---	---

Changes related to "Installation, Backup, and Upgrade" (page 27):

Object	Description	Document
System Update from the Previous Software Version (SR1075.8.3.25) (page 31)	New Requirement for the system update from the previous software version. This requirement will be updated for each new software version with the system update tasks available for the new software version.	1.0
Data Migration from the Previous Software Version (SR1075.8.3.26) (page 31)	New Requirement for the data migration from the previous software version. This requirement will be updated for each new software version with the data migration tasks available for the new software version.	1.0
Installation, Backup, and Upgrade (SR1095.4+) (page 27)	Update Renamed since the grouping requirement also covers requirements related to backup and upgrade.	1.0
Configuration of MSDS, Pictograms, and NFPA Logo (SR3146.7.1)	Deleted Display of hazard-related data is covered by the related requirements in the Functional Requirement Specification Data Management.	1.0

Changes related to "User Experience" (page 33):

Object	Description	Document
ERP - Standard Integration Scenarios (SR1095.5.1.2.15)	Deleted The extension use case is no longer supported.	1.0
Add New Fields (SR1095.5.1.2.15.1)	Deleted The extension use case is no longer supported.	1.0
Create New Custom Integration Scenario (SR1095.5.1.2.15.2)	Deleted The extension use case is no longer supported.	1.0

Object	Description	Document
Enable Consumption Message Upload (SR1095.5.1.2.15.3)	Deleted The extension use case is no longer supported.	1.0

-
-
- Rockwell Software PharmaSuite® 8.4 - Functional Requirement Specification Non-functional Requirements
-
-

C

- Compliance Requirements (SR1095.50+) • 13
 - Access control (SR1095.50.3) • 16
 - Access privileges for pre-defined set of functions (SR1095.50.3.4.2) • 20
 - Activity status (user account) (SR1095.50.3.2.6) • 19
 - Assigning users to user groups (SR1095.50.3.2.10) • 19
 - Audit trail (access privilege) (SR1095.50.3.4.1) • 20
 - Audit trail (SR1095.50.1) • 13
 - Audit trail (user group) (SR1095.50.3.3.1) • 20
 - Audit trail (user) (SR1095.50.3.2.1) • 18
 - Authentication for system login and electronic signature (SR1095.50.3.1) • 16
 - Comments (SR1095.50.2.1.1) • 15
 - Configuring electronic signatures (SR1095.50.2.2) • 15
 - Data of an electronic signature (SR1095.50.2.4) • 16
 - Default user account settings (SR1095.50.3.2.2) • 18
 - Electronic signatures (SR1095.50.2) • 13
 - Expiry date (password) (SR1095.50.3.2.9) • 19
 - Immutable login name (SR1095.50.3.2.4) • 19
 - Invalid authentication attempt (SR1095.50.3.1.2) • 17
 - Linkage to electronic record (SR1095.50.2.1.2) • 15
 - Logging unsuccessful authentication attempts (SR1095.50.3.1.3) • 17
 - Management of access privileges (SR1095.50.3.4) • 20
 - Management of user groups (SR1095.50.3.3) • 20
 - Management of users (SR1095.50.3.2) • 18
 - Multiple logins per user (SR1095.50.3.1.5) • 17
 - Number of invalid authentication attempts (SR1095.50.3.2.7) • 19
 - Number of invalid authentication attempts exceeded (SR1095.50.3.1.4) • 17
 - Password change - after login (SR1095.50.3.6.2) • 21
 - Password change - system-forced (SR1095.50.3.6.3) • 21
 - Password change - user-initiated (SR1095.50.3.6.1) • 21
 - Password change (SR1095.50.3.6) • 21
 - Password encryption (SR1095.50.3.5.1) • 21
 - Password must not be displayed (SR1095.50.3.5.2) • 21
 - Password security (SR1095.50.3.5) • 20
 - Password strength (SR1095.50.3.2.8) • 19
 - Performing electronic signatures (SR1095.50.2.3) • 15
 - Pre-defined roles and profiles (SR1095.50.3.4.3) • 20
 - Single and double electronic signatures (SR1095.50.2.1) • 14
 - Unique login name (SR1095.50.3.2.3) • 18
 - Validity timeframe (user account) (SR1095.50.3.2.5) • 19
 - Verification during authentication (SR1095.50.3.1.1) • 16
 - Viewing access privileges of users or user groups (SR1095.50.3.4.4) • 20
 - Viewing system access history (SR1095.50.3.1.6) • 17
 - Viewing users or user groups of access privileges (SR1095.50.3.4.5) • 20
- Conventions (typographical) • 2

I

- Installation, backup, and upgrade (SR1095.4+) • 27
 - Backup (SR1075.6) • 28
 - Backup and restoring (SR1075.6.3) • 28
 - Backup of a running system (SR1075.6.1) • 28
 - Basic upgrade engine (SR1075.8.1) • 29
 - Configuration parameters of Express Edition installer (SR1095.4.1.9) • 28
 - Data migration from the previous aofware version (SR1075.8.3.26) • 31
 - Existing database (SR1095.4.1.7) • 27
 - Extensibility (SR1075.8.5+) • 32
 - Initial installation (SR1095.4.1.1) • 27
 - Installation manual (SR1095.4.1.4) • 28
 - Installation of PharmaSuite (SR1095.4.1) • 27
 - Installation result (SR1095.4.1.8) • 27
 - Logging of installation (SR1095.4.1.8.1) • 27
 - Pre-install configuration (SR3146.7) • 32

Scope of Express Edition installer (SR1095.4.1.6) • 28

System update and data migration packages
(SR1075.8.3+) • 31

System update from the previous software version
(SR1075.8.3.25) • 31

System update tasks (SR1075.8.2+) • 29

Update task for Access Privilege objects
(SR1075.8.2.15) • 31

Update task for Application objects (SR1075.8.2.12) • 30

Update task for AT Definition objects (SR1075.8.2.8) • 30

Update task for ATRow objects (SR1075.8.2.4) • 30

Update task for DSX objects (generic DSX)
(SR1075.8.2.1) • 29

Update task for EAR files (SR1075.8.2.3) • 29

Update task for FSM (Flexible State Model) objects
(SR1075.8.2.7) • 30

Update task for Image objects (SR1075.8.2.6) • 30

Update task for Library objects (SR1075.8.2.5) • 30

Update task for Message objects (SR1075.8.2.2) • 29

Update task for Report Design objects (SR1075.8.2.14) • 31

Update task for UDA Definition objects (SR1075.8.2.9) • 30

Update task for Versioning Configuration objects
(SR1075.8.2.13) • 31

Update task to report customer notifications
(SR1075.8.2.11) • 30

Upgrade - Continue execution (SR1075.8.3.5) • 29

Upgrade packages (SR1075.8.5.1) • 32

Upgrade tool (SR1075.8+) • 28

Internationalization and localization (SR1095.1+) • 23

Choice of locale (SR1095.2.2) • 25

Configurable locale-specific settings (SR1095.1.3) • 25

English (United States) (SR1095.2.3.1) • 24

Internationalized output media (SR1095.2) • 23

Locale-related requirements • 24

Locales within one installation (SR1095.2.1) • 24

Provided locales (SR1095.2.3+) • 24

Skills for localization (SR1095.1.1) • 23

Supported locales (SR1095.1.2) • 24

Tool support (SR1095.2.4) • 23

L

Locale-related requirements • 24

P

Performance Requirements (SR1095.10+) • 3

Change request - Performance (SR1095.10.1+) • 5

EBR server - Performance (SR1095.10.3+) • 9

Equipment - Performance (SR1095.10.2+) • 6

Exception Dashboard - Performance (1095.10.8+) • 10

Master recipe - Performance (1095.10.4+) • 3

Performance of batch report printing 1 (1095.10.8.2) • 11

Performance of change request execution 1
(SR1095.10.1.1) • 5

Performance of change request execution 2
(SR1095.10.1.2) • 6

Performance of ETO processing 1 (SR1095.10.5.1) • 10

Performance of filter class export 1 (SR1095.10.2.5) • 7

Performance of filter class import 1 (SR1095.10.2.6) • 8

Performance of master recipe export 1 (1095.10.4.2) • 4

Performance of master recipe import 1 (1095.10.4.1) • 4

Performance of master recipe report 1 (1095.10.4.5) • 4

Performance of master recipe saving 1 (1095.10.4.3) • 4

Performance of order explosion 1 (1095.10.4.4) • 4

Performance of order processing 1 (SR1095.10.3.1) • 9

Performance of prepare for status change 1
(1095.10.4.6) • 4

Performance of printing (SR1090.4.6) • 9

Performance of recipe execution 1 (SR1095.10.7.1) • 8

Performance of review by exception 1 (1095.10.8.1) • 11

Performance of unit class export 1 (SR1095.10.2.1) • 7

Performance of unit class import 1 (SR1095.10.2.2) • 7

Performance of usage list creation 1 (SR1095.10.1.3) • 6

Printing of barcode label and dispensing report • 9

Recipe execution server - Performance (SR1095.10.7+) • 8

TOM server - Performance (SR1095.10.5+) • 10

Printing of barcode label and dispensing report • 9

S

SR1075.6 - Backup • 28

SR1075.6.1 - Backup of a running system • 28

SR1075.6.3 - Backup and restoring • 28

SR1075.8.1 - Basic upgrade engine • 29

SR1075.8.2.1 - Update task for DSX objects (generic DSX) • 29

SR1075.8.2.11 - Update task to report customer notifications • 30

SR1075.8.2.12 - Update task for Application objects • 30

SR1075.8.2.13 - Update task for Versioning Configuration objects • 31

SR1075.8.2.14 - Update task for Report Design objects • 31

SR1075.8.2.15 - Update task for Access Privilege objects • 31

SR1075.8.2.2 - Update task for Message objects • 29

SR1075.8.2.3 - Update task for EAR files • 29

SR1075.8.2.4 - Update task for ATRow objects • 30

SR1075.8.2.5 - Update task for Library objects • 30

SR1075.8.2.6 - Update task for Image objects • 30

SR1075.8.2.7 - Update task for FSM (Flexible State Model) objects • 30

SR1075.8.2.8 - Update task for AT Definition objects • 30

SR1075.8.2.9 - Update task for UDA Definition objects • 30

SR1075.8.2+ - System update tasks • 29

SR1075.8.3.25 - System update from the previous software version • 31

SR1075.8.3.26 - Data migration from the previous software version • 31

SR1075.8.3.5 - Upgrade - Continue execution • 29

SR1075.8.3+ - System update and data migration packages • 31

SR1075.8.5.1 - Upgrade packages • 32

SR1075.8.5+ - Extensibility • 32

SR1075.8+ - Upgrade tool • 28

SR1090.4.6 - Performance of printing • 9

SR1095.1.1 - Skills for localization • 23

SR1095.1.2 - Supported locales • 24

SR1095.1.3 - Configurable locale-specific settings • 25

SR1095.1+ - Internationalization and localization • 23

SR1095.10.1.1 - Performance of change request execution 1 • 5

SR1095.10.1.2 - Performance of change request execution 2 • 6

SR1095.10.1.3 - Performance of usage list creation 1 • 6

SR1095.10.1+ - Change request - Performance • 5

SR1095.10.2.1 - Performance of unit class export 1 • 7

SR1095.10.2.2 - Performance of unit class import 1 • 7

SR1095.10.2.5 - Performance of filter class export 1 • 7

SR1095.10.2.6 - Performance of filter class import 1 • 8

SR1095.10.2+ - Equipment - Performance • 6

SR1095.10.3.1 - Performance of order processing 1 • 9

SR1095.10.3+ - EBR server - Performance • 9

SR1095.10.4.1 - Performance of master recipe import 1 • 4

SR1095.10.4.2 - Performance of master recipe export • 4

SR1095.10.4.3 - Performance of master recipe saving 1 • 4

SR1095.10.4.4 - Performance of order explosion 1 • 4

SR1095.10.4.5 - Performance of master recipe report 1 • 4

SR1095.10.4.6 - Performance of prepare for status change 1 • 4

SR1095.10.4+ - Master recipe - Performance • 3

SR1095.10.5.1 - Performance of ETO processing 1 • 10

SR1095.10.5+ - TOM server - Performance • 10

SR1095.10.7.1 - Performance of recipe execution 1 • 8

SR1095.10.7+ - Recipe execution - Performance • 8

SR1095.10.8.1 - Performance of review by exception 1 • 11

SR1095.10.8.2 - Performance of batch report printing 1 • 11

SR1095.10.8+ - Exception Dashboard - Performance • 10

SR1095.10+ - Performance Requirements • 3

SR1095.2 - Internationalized output media • 23

SR1095.2.1 - Locales within one installation • 24

SR1095.2.2 - Choice of locale • 25

SR1095.2.3.1 - English (United States) • 24

SR1095.2.3+ - Provided locales • 24

SR1095.2.4 - Tool support • 23

SR1095.4.1 - Installation of PharmaSuite • 27

SR1095.4.1.1 - Initial installation • 27	SR1095.5.7.1 - GUI adaptation • 33
SR1095.4.1.4 - Installation manual • 28	SR1095.5.7.2 - Dialog boxes • 33
SR1095.4.1.6 - Scope of Express Edition installer • 28	SR1095.5.7.3 - GUI design • 33
SR1095.4.1.7 - Existing database • 27	SR1095.5.7.4 - Grid components • 34
SR1095.4.1.8 - Installation result • 27	SR1095.5+ - User experience • 33
SR1095.4.1.8.1 - Logging of installation • 27	SR1095.50.1 - Audit trail • 13
SR1095.4.1.9 - Configuration parameters of Express Edition installer • 28	SR1095.50.2 - Electronic signatures • 13
SR1095.4+ - Installation, backup, and upgrade • 27	SR1095.50.2.1 - Single and double electronic signatures • 14
SR1095.5.1.1 - Messages • 33	SR1095.50.2.1.1 - Comments • 15
SR1095.5.1.2.1 - System - Customization • 34	SR1095.50.2.1.2 - Linkage to electronic record • 15
SR1095.5.1.2.10 - System - Batch status • 35	SR1095.50.2.2 - Configuring electronic signatures • 15
SR1095.5.1.2.11 - System - ERP BOM attribute • 34	SR1095.50.2.3 - Performing electronic signatures • 15
SR1095.5.1.2.13 - Master recipe - Approval workflow • 35	SR1095.50.2.4 - Data of an electronic signature • 16
SR1095.5.1.2.14 - Expression editor - Context data • 35	SR1095.50.3 - Access control • 16
SR1095.5.1.2.16.1 - Master recipe report variants • 36	SR1095.50.3.1 - Authentication for system login and electronic signature • 16
SR1095.5.1.2.16.2 - Master workflow report variants • 36	SR1095.50.3.1.1 - Verification during authentication • 16
SR1095.5.1.2.16.3 - Batch report variants • 36	SR1095.50.3.1.2 - Invalid authentication attempt • 17
SR1095.5.1.2.16.4 - Workflow report variants • 36	SR1095.50.3.1.3 - Logging unsuccessful authentication attempts • 17
SR1095.5.1.2.16+ - Reporting • 36	SR1095.50.3.1.4 - Number of invalid authentication attempts exceeded • 17
SR1095.5.1.2.2 - Execution - Autostart • 35	SR1095.50.3.1.5 - Multiple logins per user • 17
SR1095.5.1.2.3 - System - View model-related objects • 34	SR1095.50.3.1.6 - Viewing system access history • 17
SR1095.5.1.2.4 - Building block - Name spaces • 35	SR1095.50.3.2 - Management of users • 18
SR1095.5.1.2.5 - Execution - Select order grid • 36	SR1095.50.3.2.1 - Audit trail (user) • 18
SR1095.5.1.2.6 - Building block - Mainstream GUI programming • 35	SR1095.50.3.2.10 - Assigning users to user groups • 19
SR1095.5.1.2.8 - Execution - EBR header bar • 36	SR1095.50.3.2.2 - Default user account settings • 18
SR1095.5.1.2.9 - System - Material attribute • 34	SR1095.50.3.2.3 - Unique login name • 18
SR1095.5.1.2+ - Extensibility • 34	SR1095.50.3.2.4 - Immutable login name • 19
SR1095.5.1.3.1 - Fixed point numbers • 36	SR1095.50.3.2.5 - Validity timeframe (user account) • 19
SR1095.5.1.3.2 - Floating point arithmetic • 37	SR1095.50.3.2.6 - Activity status (user account) • 19
SR1095.5.1.3.3 - Errors at calculation • 37	SR1095.50.3.2.7 - Number of invalid authentication attempts • 19
SR1095.5.1.3.4 - Input data errors • 37	SR1095.50.3.2.8 - Password strength • 19
SR1095.5.1.3.5 - Rounding • 37	SR1095.50.3.2.9 - Expiry date (password) • 19
SR1095.5.1.3+ - Precision of numerical data • 36	SR1095.50.3.3 - Management of user groups • 20
SR1095.5.1.4.2 - Touch screen • 33	SR1095.50.3.3.1 - Audit trail (user group) • 20
SR1095.5.1.4.2.1 - Virtual keyboard • 33	SR1095.50.3.4 - Management of access privileges • 20
SR1095.5.1.4+ - Device-specific usability • 33	SR1095.50.3.4.1 - Audit trail (access privilege) • 20
SR1095.5.1.5 - Concurrent data access - Pessimistic locking • 37	
SR1095.5.7 - Look and feel • 33	

SR1095.50.3.4.2 - Access privileges for pre-defined set of functions • 20
 SR1095.50.3.4.3 - Pre-defined roles and profiles • 20
 SR1095.50.3.4.4 - Viewing access privileges of users or user groups • 20
 SR1095.50.3.4.5 - Viewing users or user groups of access privileges • 20
 SR1095.50.3.5 - Password security • 20
 SR1095.50.3.5.1 - Password encryption • 21
 SR1095.50.3.5.2 - Password must not be displayed • 21
 SR1095.50.3.6 - Password change • 21
 SR1095.50.3.6.1 - Password change - user-initiated • 21
 SR1095.50.3.6.2 - Password change - after login • 21
 SR1095.50.3.6.3 - Password change - system-forced • 21
 SR1095.50+ - Compliance Requirements • 13
 SR1095.8 - Restart and recovery • 37
 SR3146.7 - Pre-install Configuration • 32

U

User experience (SR1095.5+) • 33
 Batch report variants (SR1095.5.1.2.16.3) • 36
 Building block - Mainstream GUI programming (SR1095.5.1.2.6) • 35
 Building Block - Name spaces (SR1095.5.1.2.4) • 35
 Concurrent data access - Pessimistic locking (SR1095.5.1.5) • 37
 Device-specific usability (SR1095.5.1.4+) • 33
 Dialog boxes (SR1095.5.7.2) • 33
 Errors at calculation (SR1095.5.1.3.3) • 37
 Execution - Autostart (SR1095.5.1.2.2) • 35
 Execution - EBR header bar (SR1095.5.1.2.8) • 36
 Execution - Select order grid (SR1095.5.1.2.5) • 36
 Expression editor - Context data (SR1095.5.1.2.14) • 35
 Extensibility (SR1095.5.1.2+) • 34
 Fixed point numbers (SR1095.5.1.3.1) • 36
 Floating point arithmetic (SR1095.5.1.3.2) • 37
 Grid components (SR1095.5.7.4) • 34
 GUI adaptation (SR1095.5.7.1) • 33
 GUI design (SR1095.5.7.3) • 33
 Input data errors (SR1095.5.1.3.4) • 37
 Look and feel (SR1095.5.7) • 33

Master recipe - Approval workflow (SR1095.5.1.2.13) • 35
 Master recipe report variants (SR1095.5.1.2.16.1) • 36
 Master workflow report variants (SR1095.5.1.2.16.2) • 36
 Messages (SR1095.5.1.1) • 33
 Precision of numerical data (SR1095.5.1.3+) • 36
 Reporting (SR1095.5.1.2.16+) • 36
 Restart and recovery (SR1095.8) • 37
 Rounding (SR1095.5.1.3.5) • 37
 System - Batch status (SR1095.5.1.2.10) • 35
 System - Customization (SR1095.5.1.2.1) • 34
 System - ERP BOM attribute (SR1095.5.1.2.11) • 34
 System - Material attribute (SR1095.5.1.2.9) • 34
 System - View model-related objects (SR1095.5.1.2.3) • 34
 Touch screen (SR1095.5.1.4.2) • 33
 Virtual keyboard (SR1095.5.1.4.2.1) • 33
 Workflow report variants (SR1095.5.1.2.16.4) • 36